

TrichinEasy
cod. 90060

&

PLYTricons
cod. 12070

EASY, RELIABLE, ACCURATE



Anisakis



**Trichineasy,
a complete system
for Anisakis
detection and
analysis**



Trichineasy Anisakis technical features

Operator's safety: all dedicated reagents disposable kit are predosed, in sealed containers with low concentration of HCl.

Easy to use, thanks to the ergonomic container.

Reliable results, on standardized procedures.

Repeatability granted.

High quality standards, 99% of digestion granted instead of 5% of non-digestion possibility indicated in the Italian law 2075/05.

Quick results: thanks to the original material used, the minimum quantities of liquids required and the special rotation technology of the blades digestion will be completed in **20 minutes!**

The complete system is made of:

- Digestion and filtering system;
- Analysis Kit (**Plytricons 12070**): pre dosed HCl, pre dosed pepsin 10g, cleaning solution and plastic pliers. (See photo)

Safety:

- LVD 2006/95/CE
- EN 61326:1997 + A1:1998 + A2:2001 + A3:2003
- EMC 2004/108/CE

Optional:

Fluorescent microscope.



Plytricons complete kit, 30 runs per box. [see photo]

Stereomicroscope larvae research. [see photo]

